RAPIDEC[®] ur

IVD

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RAPIDEC[®] ur

IVD

Detection of common urinary pathogens

SUMMARY AND EXPLANATION

RAPIDEC ur is a standardized micromethod for the 2-hour detection of the main pathogens responsible for urinary infections :

Escherichia coli, *Proteus mirabilis*, other *Proteus/Providencia*, *Klebsiella* group/*Enterobacter*, *Serratia marcescens*, enterococci, *Staphylococcus aureus/epidermidis*, *Staphylococcus saprophyticus* and *Candida albicans*.

PRINCIPLE

The RAPIDEC ur strip consists of 4 cupules allowing the performance of 8 tests for the detection of common urinary pathogens.

The biochemical tests are read and interpreted according to the sequence of the reactions after 2 hours of incubation at 37°C.

CONTENT OF THE KIT (Kit for 30 tests) :

- 10 strips (3 tests each)
- 10 lids
- 1 package insert

COMPOSITION OF THE STRIP

The composition of the RAPIDEC ur strip is given below in the list of tests :

CUP.	TESTS	SUBSTRATES	QTY (mg/cup.)	REACTIONS
C	Turbidity control (4 McFarland)			
S	Preparation of the bacterial suspension			
1	LANA	L-alanyl-L-alanine-4-nitroanilide	0.021	Alanine Arylamidase
2	β-GUR	4-nitrophenyl-βD-glucuronide	0.032	β-GlucURonidase
3	PPA	4-nitrophenylalanine	0.024	Phenylalanine deaminase
	PRO	2-naphthyl-βD-galactopyranoside	0.092	PROline arylamidase
	β-GAL	L-proline-7-amido-4-methylcoumarin	0.0088	β-GALactosidase
4	βG	4-methylumbelliferyl-βD-glucopyranoside	0.018	β-Glucosidase
	βX	4-methylumbelliferyl-βD-xylopyranoside	0.0156	β-Xylosidase
	IND	L-tryptophane	0.04	INDole production

- The quantities indicated may be adjusted depending on the titer of the raw materials used.
- Certain cupules contain products of animal origin, notably peptones.

REAGENTS AND MATERIAL REQUIRED BUT NOT PROVIDED**Reagents :**

- API[®] Suspension Medium, 2 ml (Ref. 70 700)
- JAMES reagent (Ref. 70 542)
- FVB reagent (Ref. 70 500)
- Oxidase (Ref. 55 635*)
- *reference not sold in certain countries : use an equivalent reagent.
- Hydrogen peroxide (3 %)

Material :

- Wooden or glass applicator sticks
- Pipettes or PSIpettes, or precision pipette
- Ampule rack
- Ampule protector
- General microbiology laboratory equipment including Ultra violet lamp (365 nm)

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline* - Current revision". For additional handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories - CDC/NIH - Latest edition", or to the regulations currently in use in each country.
- Do not use reagents past the expiration date.
- Before use, check that the packaging of the various components is intact.
- Do not use strips which have been damaged : cupules deformed, desiccant sachet open, etc.
- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.
- Interpretation of the test results should be made taking into consideration the patient history, the source of the specimen, colonial and microscopic morphology of the strain and, if necessary, the results of any other tests performed.

STORAGE CONDITIONS

The strips should be stored at 2-8°C until the expiration date indicated on the packaging.

SPECIMENS (COLLECTION AND PREPARATION)

RAPIDEC ur is not for use directly with clinical specimens. The microorganisms to be identified must first be isolated on a suitable culture medium according to standard microbiological techniques.

INSTRUCTIONS FOR USE

Preparation of the strip

- Remove the strip from its packaging.
- Write the strain or specimen reference numbers on the strip next to the letters A, B or C.
- Place the lid on the strip.
If the strip is not completely used, cut off the number of tests required using a pair of scissors. Replace the remainder of the strip in the refrigerator in its packaging and use it within 2 days.

Preparation of the inoculum

- Open an ampule of API® Suspension Medium (2 ml) as indicated in the paragraph "Warnings and Precautions" of the package insert for this product, or use any tube containing 2 ml of distilled water without additives : dispense 250 µl into each of the cupules C and S. The level of the water must be 1 mm below the top edge of the cupule.
- Gently homogenize the contents of cupule C with an applicator stick. Proceed immediately to the next step.
- Using the end of a new stick, pick up several colonies of the same morphology and homogenize them in cupule S until a turbidity equivalent to that in cupule C (4 McFarland) is obtained, i.e. about 2 to 5 colonies. It is recommended to use young cultures (18-24 hours old).

NOTES :

- Do not use saline medium for rehydration of cupules C and S (it would prevent the appearance of turbidity in cupule C).
- Any spatter of liquid from cupule C into cupules S, 1, 2, 3 or 4 will kill the bacteria and prevent the enzymatic reactions.
- It is easier to compare turbidities by placing the strip on a black background.

Inoculation of the strip

- Transfer 50 µl of the bacterial suspension (S) into the corresponding cupules 1, 2, 3 and 4.
- Put the excess suspension back into cupule S.
- Withdraw the liquid from cupule C to avoid contamination of the other cupules.
- Place a lid on the strip.
- Incubate the RAPIDEC ur strip for 2 - 2 ¼ hours at 36°C ± 2°C.

READING AND INTERPRETATION

Reading the strip

Refer to the Reading Table.

- Only the spontaneous colored reactions should be read systematically (LANA, β-GUR and PPA).
- The reactions requiring a UV lamp (PRO and βG + βX) and reagents (IND, β-GAL, CAT et OX) should only be read if indicated in the Reading Table.
- Reactions requiring a reagent :
 - INDole test (IND)
Add 1 drop of JAMES reagent to cupule 4. The reaction takes place immediately.
 - β-GALactosidase test (β-GAL)
Add 1 drop of FVB reagent to cupule 3. Wait 5 minutes before reading the reaction.
 - CATalase test (CAT)
Add 1 drop of hydrogen peroxide (3 %) to cupule 2. Read the reaction after 1-2 minutes.
 - OXidase test (OX)
The oxidase test must be performed according to the manufacturer's instructions for use.

NOTES :

- The incubation time can exceed 2 hours, but should not exceed 4 hours. After this length of time, it is difficult to interpret the reactions and the results could be false.
- Reading is made easier by placing the strip on a white background.

Interpretation

Refer to the Reading Table.

Continuation of the analysis

After reading, the bacterial suspension which remains in cupule S can be used for complete identification (API® strip), an antimicrobial susceptibility test (ATB® strip) or a subculture.

QUALITY CONTROL

The media, strips and reagents are systematically controlled at various stages of their manufacture. For those users who wish to perform their own quality control tests with the strip, it is preferable to use the strain **1. *Escherichia coli* ATCC® 10536** or else one of the following strains :

- | | | | |
|--------------------------------|-------------|--|-------------|
| 2. <i>Providencia stuartii</i> | ATCC® 25825 | 4. <i>Staphylococcus saprophyticus</i> | ATCC® 43867 |
| 3. <i>Serratia marcescens</i> | ATCC® 43862 | | |

ATCC : American Type Culture Collection, 10801 University Boulevard, Manassas, VA 20110-2209, USA.

	Spontaneous reactions			Fluorescent reactions		Reactions requiring a reagent		
	Cup 1	Cup 2	Cup 3	Cup 3	Cup 4	Cup 2	Cup 3	Cup 4
	LANA	β-GUR	PPA	PRO	βG + βX	CAT	β-GAL	IND
1.	+	+	–	– (1)	– (1)		+	+
2.	+	–	+	– (1)	– (1)		– (1)	+
3.	+	–	–	+	+		+	– (1)
4.	–	–	–	–	–	+	+	

(1) These reactions are not necessary for the detection of these species (see Reading Table)

It is the responsibility of the user to perform Quality Control in accordance with any local applicable regulations.

LIMITATIONS OF THE METHOD

- RAPIDEC ur is designed to detect bacteria which cause urinary infections. It is not considered as an identification system. In case of doubt, complete identification should be performed.
- The following isolation media can be used : blood agar, MacConkey, BCP, CLED, Hektoen and DCLS media, and, to a lesser degree, SS Medium, prevent the IND reaction and therefore must not be used.

- Certain isolation media cause a pigmentation of the colonies. This does generally not affect the reading of reactions and these media can therefore be used. The EMB medium, however, contains an indicator which can mask certain reactions and make it difficult to read the results. This must be taken into account when performing the reading.
- Only pure cultures of a single organism should be used.

RANGE OF EXPECTED RESULTS

Consult the Reading Table at the end of this package insert for the range of expected results for the various biochemical reactions.

PERFORMANCE

Species	Number of strains	Taxon		
		Correct	G(-) OX (-) rod	Incorrect
<i>Escherichia coli</i>	73	72	1	0
<i>Klebsiella pneumoniae</i>	37	37	0	0
<i>Proteus mirabilis</i>	37	37	0	0
<i>Pseudomonas aeruginosa</i>	15	15	0	0
<i>Enterobacter cloacae</i>	4	4	0	0
<i>Proteus vulgaris</i>	4	4	0	0
<i>Citrobacter freundii</i>	3	N/A	3	0
<i>Acinetobacter</i> spp	2	N/A	2	0
<i>Morganella morganii</i>	1	1	0	0
<i>Providencia stuartii</i>	1	1	0	0
<i>Citrobacter diversus</i>	1	1	0	0
<i>Klebsiella oxytoca</i>	1	1	0	0
<i>Staphylococcus aureus</i> & <i>S. epidermidis</i>	19	19	0	0
<i>Staphylococcus saprophyticus</i>	7	70	0	0
<i>Enterococcus</i> spp	24	24	0	0
<i>Candida albicans</i>	8	8	0	0
TOTAL	237	231	6	0
%		97.5	2.5	0.0

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their type and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

READING TABLE

SPONTANEOUS REACTIONS			OTHER REACTIONS				RESULTS
① LANA	② β-GUR	③ PPA					
yellow	yellow	colorless or beige					<i>E. coli</i>
yellow	colorless	orange or brown	④ JAMES		→ colorless		<i>P. mirabilis</i>
			IND		→ pink		<i>Proteus vulgaris</i> <i>Morganella</i> <i>Providencia</i>
yellow	colorless	colorless or beige	③ UV PRO	④ UV βG-βX			
			fluorescent	fluorescent	OX Test (1)	→ positive	<i>Pseudomonas</i>
					OX Test (1)	→ negative	<i>Serratia</i>
			colorless	fluorescent	④ JAMES IND	→ pink	<i>K. oxytoca</i> <i>C. diversus</i> (2)
						→ colorless	<i>K. pneumoniae</i> <i>Enterobacter</i> <i>Citrobacter</i> (2)
			colorless or fluorescent	colorless	④ JAMES IND	→ pink	<i>E. coli</i> / β-GUR (-)
						→ colorless :	
						OX Test (1) → positive	<i>Pseudomonas</i>
						OX Test (1) → negative	Gram (-) OX (-) rod
colorless or very pale yellow	colorless	colorless or beige	③ UV PRO	④ UV βG-βX	② H ₂ O ₂ CAT		
			colorless	fluorescent	no bubbles		
			colorless	colorless	bubbles	③ FVB β-GAL	→ purple <i>S. saprophyticus</i>
							→ pinkish beige <i>S. aureus</i> <i>S. epidermidis</i>
			fluorescent	colorless	bubbles		
							<i>C. albicans</i>

① indicates the number of the cupule.

(1) Perform the test according to the manufacturer's instructions for use.

(2) The motility test can easily separate the genus *Klebsiella* (-) from the genera *Enterobacter* (+) and *Citrobacter* (+).

To carry out this test :

- add 200 µl of Brain-Heart broth to cupule S.

- incubate for approximately 2 hours at 36°C ± 2°C.

- observe the motility under a microscope using the hanging drop technique.

If the combination of colors obtained is not listed in the table, it may be due to a mixture of organisms or to species not normally encountered in urinary tract infections. Check the purity of the strain and inoculate an API® identification strip.

PROCEDURE
LITERATURE REFERENCES
INDEX OF SYMBOLS

p. I
p. II
p. III

bioMérieux® SA

au capital de 11 879 045 €
673 620 399 RCS LYON
69280 Marcy-l'Etoile / France
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>

bioMérieux, Inc

Box 15969,
Durham, NC 27704-0969 / USA
Tel. (1) 919 620 20 00
Fax (1) 919 620 22 11

Printed in France

The logo is a registered and protected trademark of bioMérieux SA or one of its subsidiaries.

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